

Product Name: GAD67 Rabbit Polyclonal Antibody
Catalog #: APRab00173



Summary

Production Name	GAD67 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	GAD1
Alternative Names	Glutamate decarboxylase 1; 67 kDa glutamic acid decarboxylase; GAD 67; GAD1; GAD; GAD67
Gene ID	2571
SwissProt ID	Q99259.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 67 kDa; Observed MW: 67 kDa

Background

Product Name: GAD67 Rabbit Polyclonal Antibody
Catalog #: APRab00173




Catalyzes the production of GABA.

Research Area

Signal Transduction

Image Data



Western blot analysis of GAD67 in (1) Jurkat lysates; (2) C6 lysates using GAD67 antibody.

Note

For research use only.