

Product Name: LMO2 Rabbit Polyclonal Antibody**Catalog #: APRab00208**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,FC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 18 kDa; Observed MW: 18 kDa

Antigen Information

Gene Name	LMO2
Alternative Names	TTG2; LMO-2; RBTN2; RHOM2; RBTNL1
Gene ID	4005
SwissProt ID	P25791
Immunogen	Recombinant protein of human LMO2


Background

Acts with TAL1/SCL to regulate red blood cell development. Also acts with LDB1 to maintain erythroid precursors in an immature state.

Research Area

Cardiovascular

Image Data



Western blot analysis of LMO2 in Raji lysates using LMO2 antibody.