

Product Name: TCTP Rabbit Polyclonal Antibody**Catalog #: APRab00222**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Antigen Information

Gene Name	TPT1
Alternative Names	Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23
Gene ID	7178
SwissProt ID	P13693
Immunogen	A synthetic peptide of human TCTP


Background

Involved in calcium binding and microtubule stabilization.

Research Area

Immunology

Image Data



Western blot analysis of TCTP Receptor in HeLa lysates using TCTP antibody.