

Product Name: Follistatin Rabbit Polyclonal Antibody**Catalog #: APRab00231**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Antigen Information

Gene Name	FST
Alternative Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene
Gene ID	10468
SwissProt ID	P19883
Immunogen	Recombinant protein of human Follistatin


Background

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Research Area

Signal Transduction

Image Data



Western blot analysis of Follistatin in HepG2 lysates using Follistatin antibody.