

Product Name: Factor XIIIa Rabbit Polyclonal Antibody
Catalog #: APRab00235



Summary

Production Name	Factor XIIIa Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	F13A1
Alternative Names	Coagulation factor XIIIa; F13A; F13a1; Fibrin stabilizing factor; A subunit; Fibrinolygase; TGase; Transglutaminase A chain
Gene ID	2162
SwissProt ID	P00488.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 FC: 1:50-1:100
Molecular Weight	Calculated MW: 83 kDa; Observed MW: 83 kDa

Background

Product Name: Factor XIIIa Rabbit Polyclonal Antibody
Catalog #: AP Rab00235




Factor XIII is activated by thrombin and calcium ion to a transglutaminase that catalyzes the formation of gamma-glutamyl-epsilon-lysine cross-links between fibrin chains, thus stabilizing the fibrin clot. Also cross-link alpha-2-plasmin inhibitor, or fibronectin, to the alpha chains of fibrin.

Research Area

Cardiovascular

Image Data



Western blot analysis of Factor XIIIa in Human spleen lysates using Factor XIIIa antibody.

Note

For research use only.