

Product Name: StAR Rabbit Polyclonal Antibody
Catalog #: APRab00244



Summary

Production Name	StAR Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	STAR
Alternative Names	STARD1
Gene ID	6770
SwissProt ID	P49675.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 32 kDa; Observed MW: 32 kDa

Background

Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer

Product Name: StAR Rabbit Polyclonal Antibody
Catalog #: APRab00244




membrane to the inner mitochondrial membrane.

Research Area

Signal Transduction

Image Data



Western blot analysis of StAR in SKBR3 lysates using StAR antibody.

Note

For research use only.