

---

**Product Name: SFPQ Rabbit Polyclonal Antibody****Catalog #: APRab00314**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification</b>	Affinity Chromatography

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
<b>Molecular Weight</b>	Calculated MW: 76 kDa; Observed MW: 100 kDa

**Antigen Information**

<b>Gene Name</b>	SFPQ
<b>Alternative Names</b>	PSF; POMP100; PPP1R140
<b>Gene ID</b>	6421
<b>SwissProt ID</b>	P23246
<b>Immunogen</b>	A synthetic peptide of human SFPQ

**Background**


DNA- and RNA binding protein, involved in several nuclear processes. Essential pre-mRNA splicing factor required early in spliceosome formation and for splicing catalytic step II, probably as a heteromer with NONO. Binds to pre-mRNA in

spliceosome C complex, and specifically binds to intronic polypyrimidine tracts.

## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of SFPQ in A431 lysates using SFPQ antibody.