

Product Name: RBPJK Rabbit Polyclonal Antibody
Catalog #: APRab00318



Summary

Production Name	RBPJK Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	RBPJ
Alternative Names	SUH; csl; AOS3; CBF1; KBF2; RBP-J; RBPJK; IGKJRB; RBPSUH; IGKJRB1
Gene ID	3516
SwissProt ID	Q06330.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 60 kDa

Background

Transcriptional regulator that plays a central role in Notch signaling, a signaling pathway involved in cell-cell

Product Name: RBPJK Rabbit Polyclonal Antibody
Catalog #: APRab00318




communication that regulates a broad spectrum of cell-fate determinations. Acts as a transcriptional repressor when it is not associated with Notch proteins.

Research Area

Neuroscience

Image Data



Western blot analysis of RBPJK in MCF-7 lysates using RBPJK antibody.

Note

For research use only.