
Product Name: GST3 Rabbit Polyclonal Antibody**Catalog #: APRab00343**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Antigen Information

Gene Name	GSTP1
Alternative Names	GSTP1; FAEES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1
Gene ID	2950
SwissProt ID	P09211
Immunogen	A synthetic peptide of human GST3/GST pi


Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Research Area

Cardiovascular

Image Data



Western blot analysis of GST3 / GST pi in PC-3 lysates using GST3 antibody.