

Product Name: Fas Ligand Rabbit Polyclonal Antibody**Catalog #: APRab00379**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 31 kDa; Observed MW: 31 kDa

Antigen Information

Gene Name	FASLG FASLG; APT1LG1; CD95L; FASL; TNFSF6; Tumor necrosis factor ligand superfamily member 6;
Alternative Names	Apoptosis antigen ligand; APTL; CD95 ligand; CD95-L; Fas antigen ligand; Fas ligand; FasL; CD antigen CD178
Gene ID	356
SwissProt ID	P48023
Immunogen	


Background

Cytokine that binds to TNFRSF6/FAS, a receptor that transduces the apoptotic signal into cells.

Research Area

Cell Biology

Image Data



Western blot analysis of Fas Ligand in various lysates using Fas Ligand antibody.