
Product Name: AMBN Rabbit Polyclonal Antibody**Catalog #: APRab00425**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 48 kDa; Observed MW: 48 kDa

Antigen Information

Gene Name	AMBN
Alternative Names	A11F; Ameloblastin
Gene ID	258
SwissProt ID	Q9NP70
Immunogen	


Background

Involved in the mineralization and structural organization of enamel.

Research Area

Signal Transduction

Image Data



Western blot analysis of AMBN in HeLa lysates using Ameloblastin antibody.