

Product Name: Lactate Dehydrogenase D Rabbit Polyclonal Antibody
Catalog #: APRab00438

Summary

Production Name	Lactate Dehydrogenase D Rabbit Polyclonal Antibody
Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	LDHD
Alternative Names	Probable D-lactate dehydrogenase; mitochondrial; DLD; Lactate dehydrogenase D
Gene ID	197257.0
SwissProt ID	Q86WU2.

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 55 kDa; Observed MW: 55 kDa

**Product Name: Lactate Dehydrogenase D Rabbit
Polyclonal Antibody**
Catalog #: APRab00438


Background

Involved in D-lactate, but not L-lactate catabolic process.

Research Area

Signal Transduction

Image Data



Western blot analysis of Lactate Dehydrogenase D in Hela, mouse Skeletal Muscle, mouse Liver, rat Kidney lysates using LDHD antibody.

Note

For research use only.