
Product Name: PEPC Rabbit Polyclonal Antibody**Catalog #: APRab00439**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Plants
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 105-109 kDa

Antigen Information

Gene Name	PGC
Alternative Names	PGC; Gastricsin; Pepsinogen C
Gene ID	5225
SwissProt ID	P20142
Immunogen	Recombinant protein of human Pepsinogen II/PGC


Background

Hydrolyzes a variety of proteins.

Research Area

Cardiovascular

Image Data



Western blot analysis of PEPC in Maize lysates using PEPC antibody.