

Product Name: Bestrophin 2 Rabbit Polyclonal Antibody**Catalog #: APRab00442**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 57 kDa; Observed MW: 50 kDa

Antigen Information

Gene Name	BEST2
Alternative Names	BEST 2 antibody; FLJ20132 antibody; Vitelliform macular dystrophy 2 homolog antibody
Gene ID	54831
SwissProt ID	Q8NFU1
Immunogen	


Background

Forms calcium-sensitive chloride channels.

Research Area

Signal Transduction

Image Data



Western blot analysis of Bestrophin 2 in mouse Brain, rat Brain lysates using Bestrophin2 antibody.