

Product Name: Neuromedin B Receptor Rabbit Polyclonal Antibody**Catalog #: APRab00446**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 50 kDa

Antigen Information

Gene Name	NMBR
Alternative Names	Neuromedin B receptor antibody; Bombesin Receptor 1 antibody; BB1 antibody
Gene ID	4829
SwissProt ID	P28336
Immunogen	


Background

Receptor for neuromedin-B.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Neuromedin B Receptor in Human Brain, rat Brain lysates using NMBR(Bombesin Receptor)antibody.