

Product Name: Epigen Rabbit Polyclonal Antibody**Catalog #: APRab00572**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 17 kDa; Observed MW: 17 kDa

Antigen Information

Gene Name	EPGN
Alternative Names	EPGN; Epigen; Epithelial mitogen; EPG
Gene ID	255324
SwissProt ID	Q6UW88
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human EPGN. AA range:11-60


Background

Promotes the growth of epithelial cells.


Research Area

Signal Transduction

Image Data



Western blot analysis of Epigen in 293T, 823 lysates using Epigen antibody.



Western blot analysis of Epigen in 293T lysates using Epigen antibody.