
Product Name: SOX17 Rabbit Polyclonal Antibody**Catalog #: APRab00599**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 44 kDa

Antigen Information

Gene Name	SOX17
Alternative Names	SOX17; Transcription factor SOX-17
Gene ID	64321
SwissProt ID	Q9H6I2
Immunogen	Synthesized peptide derived from Sox-17 . at AA range: 90-170


Background

Acts as transcription regulator that binds target promoter DNA and bends the DNA.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of SOX17 in Jurkat lysates using Sox17 antibody.