

Product Name: Interferon epsilon Rabbit Polyclonal Antibody
Catalog #: APRab00616

Summary

Production Name	Interferon epsilon Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	IFNE
Alternative Names	IFNE; IFNE1; Interferon epsilon; IFN-epsilon; Interferon epsilon-1
Gene ID	338376
SwissProt ID	Q86WN2.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 24 kDa; Observed MW: 24 kDa

Background

Product Name: Interferon epsilon Rabbit Polyclonal Antibody
Catalog #: APRab00616




Type I interferon required for maintaining basal levels of IFN-regulated genes, including 2'-5'-oligoadenylate synthetase, IRF7 and ISG15, in the female reproductive tract.

Research Area

Immunology

Image Data



Western blot analysis of Interferon epsilon in Jurkat lysates using Interferon epsilon antibody.

Note

For research use only.