
Product Name: FGF11 Rabbit Polyclonal Antibody**Catalog #: APRab00628**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Antigen Information

Gene Name	FGF11
Alternative Names	FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3
Gene ID	2256
SwissProt ID	Q92914
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human FGF11. AA range:21-70


Background

Probably involved in nervous system development and function.

Research Area

Cardiovascular

Image Data



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.