

Product Name: TCEAL6 Rabbit Polyclonal Antibody**Catalog #: APRab00636**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 22 kDa; Observed MW: 32 kDa

Antigen Information

Gene Name	TCEAL6
Alternative Names	TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6
Gene ID	158931.0
SwissProt ID	Q6IPX3
Immunogen	


Background

May be involved in transcriptional regulation.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TCEAL6 in CoLo lysates using TCEAL6 antibody.