
Product Name: Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit Polyclonal Antibody
Catalog #: APRab00711

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,FC 1:50-1:100
Molecular Weight	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa

Antigen Information

Gene Name	PAK4/5/6
Alternative Names	PAK4; PAK5; PAK6
Gene ID	10298/57144/56924
SwissProt ID	O96013/Q9P286/Q9NQUS
Immunogen	A synthetic Phosphorylated peptide corresponding to residues target protein


Background

Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, growth, proliferation or cell survival.

Research Area

Signal Transduction

Image Data



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) in U87-MG lysates using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) antibody.