

Product Name: KCTD7 Rabbit Polyclonal Antibody**Catalog #: APRab00800**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 33 kDa; Observed MW: 33 kDa

Antigen Information


Gene Name	KCTD7
Alternative Names	BTB/POZ domain-containing protein KCTD7
Gene ID	154881
SwissProt ID	Q96MP8
Immunogen	Synthesized peptide derived from human KCTD7. at AA range: 181-230

Background

May be involved in the control of excitability of cortical neurons.

Research Area

Image Data



Western blot analysis of KCTD7 in mouse brain , mouse spleen lysates using KCTD7 antibody.