

Product Name: TNFAIP8L3 Rabbit Polyclonal Antibody
Catalog #: APRab00801



Summary

Production Name	TNFAIP8L3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	TNFAIP8L3
Alternative Names	Tumor necrosis factor alpha-induced protein 8-like protein 3; TNF alpha-induced protein 8-like protein 3; TNFAIP8-like protein 3
Gene ID	388121
SwissProt ID	Q5GJ75.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 33 kDa; Observed MW: 33 kDa

Background

Product Name: TNFAIP8L3 Rabbit Polyclonal Antibody
Catalog #: APRab00801




Acts as a lipid transfer protein.

Research Area

Immunology

Image Data



Western blot analysis of TNFAIP8L3 in 293T , mouse brain lysates using TNFAIP8L3 antibody.

Note

For research use only.