

Product Name: Collagen XIX alpha 1 Rabbit Polyclonal Antibody
Catalog #: APRab00803



Summary

Production Name	Collagen XIX alpha 1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	COL19A1
Alternative Names	Collagen alpha-1(XIX) chain (Collagen alpha-1(Y) chain)
Gene ID	1310
SwissProt ID	Q14993.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 115 kDa; Observed MW: 130 kDa

Background


Product Name: Collagen XIX alpha 1 Rabbit Polyclonal Antibody
Catalog #: APRab00803

May act as a cross-bridge between fibrils and other extracellular matrix molecules.

Research Area

Cell Biology

Image Data



Western blot analysis of Collagen XIX alpha 1 in HELA lysates using Collagen XIX alpha 1 antibody.

Note

For research use only.