

Product Name: C8 alpha Rabbit Polyclonal Antibody
Catalog #: APRab00807



Summary

Production Name	C8 alpha Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	C8A
Alternative Names	Complement component C8 alpha chain (Complement component 8 subunit alpha)
Gene ID	731
SwissProt ID	P07357.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 65 kDa; Observed MW: 60 kDa

Background

Constituent of the membrane attack complex (MAC) that plays a key role in the innate and adaptive immune response by

Product Name: C8 alpha Rabbit Polyclonal Antibody
Catalog #: APRab00807




forming pores in the plasma membrane of target cells.

Research Area

Immunology

Image Data



Western blot analysis of C8 alpha in mouse kidney, mouse heart lysates using C8 alpha antibody.

Note

For research use only.