

Product Name: Acetyl-Lysine Rabbit Polyclonal Antibody**Catalog #: APRab00848**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Acetylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Refer to figures

Antigen Information

Gene Name	-
Alternative Names	-
Gene ID	-
SwissProt ID	-
Immunogen	


Background

Acetylation of lysine residues within proteins has emerged as an important mechanism used by cells to overcome this repression.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Acetyl-Lysine in Hela, 3T3, rat Brain lysates using Acetyl Lysine antibody.