

Product Name: Human lambda Light Chain Rabbit Polyclonal Antibody
Catalog #: APRab00905

Summary

Production Name	Human lambda Light Chain Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	IGLC1
Alternative Names	IGLC1
Gene ID	3537.0
SwissProt ID	P0CG04.

Application

Dilution Ratio	WB: 1:500-1:1000 ELISA: 1:10000
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 23 kDa

Product Name: Human lambda Light Chain Rabbit Polyclonal Antibody
Catalog #: APRab00905




Background

Constant region of immunoglobulin light chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes. In the recognition phase of humoral immunity, the membrane-bound immunoglobulins serve as receptors which, upon binding of a specific antigen, trigger the clonal expansion and differentiation of B lymphocytes into immunoglobulins-secreting plasma cells.


Research Area

Immunology

Image Data



Indirect ELISA assay for Rabbit Human kappa light chain antibody. Antigen coating concentration: 2ug/ml.



Western blot analysis of Human lambda Light Chain in total Human IgG reduced and denatured in 2X SDS loading buffer with DTT using Human lambda Light Chain antibody

Note

For research use only.