

---

**Product Name: Phospho-HDAC5 (Ser498) Rabbit Polyclonal Antibody****Catalog #: APRab00916**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC
<b>Reactivity</b>	Human,Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Phosphorylated
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100
<b>Molecular Weight</b>	Calculated MW: 122 kDa; Observed MW: 122 kDa

**Antigen Information**

<b>Gene Name</b>	HDAC5
<b>Alternative Names</b>	HD5; NY-CO-9
<b>Gene ID</b>	10014
<b>SwissProt ID</b>	Q9UQL6
<b>Immunogen</b>	A synthetic Phosphorylated peptide corresponding to residues target protein

**Background**


Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle

progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes. Involved in muscle maturation by repressing transcription of myocyte enhancer MEF2C. During muscle differentiation, it shuttles into the cytoplasm, allowing the expression of myocyte enhancer factors. Involved in the MTA1-mediated epigenetic regulation of ESR1 expression in breast cancer.


## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of Phospho-HDAC5 (Ser498) in 293 lysates using Phospho-HDAC5 (Ser498) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using HDAC5(Phospho-Ser498) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.