

---

**Product Name: Phospho-OXSR1 (Thr185) Rabbit Polyclonal Antibody****Catalog #: APRab00923**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Phosphorylated
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100
<b>Molecular Weight</b>	Calculated MW: 58 kDa; Observed MW: 65 kDa

**Antigen Information**

<b>Gene Name</b>	OXSR1
<b>Alternative Names</b>	OSR1
<b>Gene ID</b>	9943
<b>SwissProt ID</b>	O95747
<b>Immunogen</b>	A synthetic Phosphorylated peptide corresponding to residues target protein


**Background**

Phosphorylates RELL1, RELL2 and RELT.


## Research Area

Signal Transduction

## Image Data



Western blot analysis of Phospho-OXR1 (Thr185) in HepG2 lysates using Phospho-OXR1 (Thr185) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using OSR1 (Phospho-Thr185) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.