
Product Name: Phospho-FGFR4 (Tyr642) Rabbit Polyclonal Antibody**Catalog #: APRab00928**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 88 kDa

Antigen Information

Gene Name	FGFR4
Alternative Names	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD334
Gene ID	2264
SwissProt ID	P22455
Immunogen	A synthetic Phosphorylated peptide corresponding to residues target protein


Background

Receptor for acidic fibroblast growth factor. Does not bind to basic fibroblast growth factor. Binds FGF19.

Research Area

Cardiovascular

Image Data



Western blot analysis of Phospho-FGFR4 (Tyr642) in Jurkat lysates using Phospho-FGFR4 (Tyr642) antibody.