
Product Name: Annexin V Rabbit Polyclonal Antibody (HRP)**Catalog #: APRab01310**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Monkey
Conjugation	HRP
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 36 kDa; Observed MW: 30 kDa

Antigen Information

Gene Name	ANXA5
Alternative Names	PP4; ANX5; ENX2; RPRGL3.
Gene ID	308
SwissProt ID	P08758
Immunogen	A synthetic peptide of human Annexin V/ANXA5


Background

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

Research Area

Cell Biology

Image Data



Western blot analysis of Annexin V in HeLa and COS7 lysates using Annexin V antibody.conjugated with HRP