

Product Name: Cytokeratin 9 Rabbit Polyclonal Antibody
Catalog #: APRab01336

Summary

Production Name	Cytokeratin 9 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	KRT9
Alternative Names	K9; CK-9; EPPK
Gene ID	3857
SwissProt ID	P35527.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 62 kDa

Background

May serve an important special function either in the mature palmar and plantar skin tissue or in the morphogenetic


Product Name: Cytokeratin 9 Rabbit Polyclonal Antibody
Catalog #: APRab01336

program of the formation of these tissues. Plays a role in keratin filament assembly.

Research Area

Cell Biology

Image Data



Western blot analysis of Cytokeratin 9 in HeLa lysates using Cytokeratin 9 antibody.

Note

For research use only.