

---

**Product Name: TDP 43 Rabbit Polyclonal Antibody****Catalog #: APRab01342**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,FC,IP
<b>Reactivity</b>	Human,Mouse,Rat,Zebrafish
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification</b>	Affinity Chromatography

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 45 kDa; Observed MW: 45 kDa

**Antigen Information**

<b>Gene Name</b>	TARDBP
<b>Alternative Names</b>	ALS10; TDP-43
<b>Gene ID</b>	23435
<b>SwissProt ID</b>	Q13148
<b>Immunogen</b>	A synthetic peptide of human TDP43

**Background**


TDP43 (TAR DNA-binding protein 43) is involved in transcriptional regulation and exon splicing. While normal TDP43 is a nuclear protein, pathological TDP43 is a component of insoluble aggregates in patients with frontotemporal lobar

degeneration (FTLD) and amyotrophic lateral sclerosis (ALS).

## Research Area

Microbiology

## Image Data



Western blot analysis of TDP43 in K562 lysates using TDP 43 antibody.