

Product Name: IKZF3 Rabbit Polyclonal Antibody
Catalog #: APRab01382



Summary

Production Name	IKZF3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	IKZF3
Alternative Names	IKZF3; ZNFN1A3; Zinc finger protein Aiolos; Ikaros family zinc finger protein 3
Gene ID	22806
SwissProt ID	Q9UKT9.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 FC: 1:50-1:100
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 70 kDa

Background

Product Name: IKZF3 Rabbit Polyclonal Antibody
Catalog #: APRab01382




Transcription factor that plays an important role in the regulation of lymphocyte differentiation. Plays an essential role in regulation of B-cell differentiation, proliferation and maturation to an effector state. Involved in regulating BCL2 expression and controlling apoptosis in T-cells in an IL2-dependent manner.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of IKZF3 in Raji, Ramos lysates using IKZF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IKZF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.