

Product Name: Human IgG Rabbit Polyclonal Antibody
Catalog #: APRab03306



Summary

Production Name	Human IgG Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	IGHG1
Alternative Names	IGHG1; Ig gamma-1 chain C region
Gene ID	-
SwissProt ID	-.

Application

Dilution Ratio	ELISA: 1:10000
Molecular Weight	-

Background

IgG is a monomeric immunoglobulin, built of two heavy chains gamma and two light chains. Each molecule has two

Product Name: Human IgG Rabbit Polyclonal Antibody
Catalog #: APRab03306




antigen binding sites. This is the most abundant immunoglobulin and is approximately equally distributed in blood and in tissue liquids, constituting 75% of serum immunoglobulins in humans.

Research Area

Immunology

Image Data



Indirect ELISA assay for antiHuman IgG rabbit pAb. Antigen coating concentration: 4ug/ml.

Note

For research use only.