

Product Name: AIFM2/FSP1 Rabbit Polyclonal Antibody
Catalog #: APRab03391



Summary

Production Name	AIFM2/FSP1 Rabbit Polyclonal Antibody
Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	AIFM2
Alternative Names	AIFM2; AMID; PRG3; Apoptosis-inducing factor 2; Apoptosis-inducing factor homologous mitochondrion-associated inducer of death; Apoptosis-inducing factor-like mitochondrion-associated inducer of death; p53-responsive gene 3 protein
Gene ID	84883
SwissProt ID	Q9BRQ8.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 38 kDa

Product Name: AIFM2/FSP1 Rabbit Polyclonal Antibody
Catalog #: APRab03391




Background

A NAD(P)H-dependent oxidoreductase involved in cellular oxidative stress response.

Research Area

Cell Biology

Image Data



Western blot analysis of FSP1 in COS7 lysates using FSP1 antibody.

Note

For research use only.