
Product Name: GABARAP Rabbit Polyclonal Antibody**Catalog #: APRab03631**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 14 kDa; Observed MW: 14,16 kDa

Antigen Information

Gene Name	GABARAP
Alternative Names	GABARAP; GABARAPL1; GABARAPL2
Gene ID	11337
SwissProt ID	O95166
Immunogen	A synthetic peptide of human GABARAP


Background

May play a role in intracellular transport of GABA(A) receptors and its interaction with the cytoskeleton.


Research Area

Neuroscience

Image Data



Western blot analysis of GABARAP in A549, 3T3 and HeLa lysates using GABARAP antibody.



Western blot analysis of GABARAP in rat Hippocampus, rat Cerebral Cortex, rat Brain and A549 lysates using GABARAP antibody.