
Product Name: SPHK1 (Phospho-Ser225) Rabbit Polyclonal Antibody**Catalog #: APRab06084**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio WB 1:500-1:2000

Molecular Weight

Antigen Information

Gene Name	SPHK1
Alternative Names	Sphingosine kinase 1 (SK 1) (SPK 1) (EC 2.7.1.91)
Gene ID	8877.0
SwissProt ID	Q9NYA1
Immunogen	Synthesized peptide derived from human SPHK1 (Phospho-Ser225)


Background

catalytic activity:ATP + sphinganine = ADP + sphinganine 1-phosphate.,catalytic activity:ATP + sphingosine = ADP + sphingosine 1-phosphate.,cofactor:Magnesium.,function:Catalyzes the phosphorylation of sphingosine to form sphingosine 1-

phosphate (SPP), a lipid mediator with both intra-and extracellular functions. Also acts on D-erythro-sphingosine and to a lesser extent sphinganine, but not other lipids, such as D,L-threo-dihydrosphingosine, N,N-dimethylsphingosine, diacylglycerol, ceramide, or phosphatidylinositol.,similarity:Contains 1 DAGKc domain.,subunit:Interacts with ACY1 (By similarity). Binds to calmodulin. Interacts with SPHKAP.,tissue specificity:Widely expressed with highest levels in adult liver, kidney, heart and skeletal muscle.,

Research Area

Image Data



Western Blot analysis of 1 MCF-7 cell, 2 LPS 100ng/mL 30min treated ,using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000