

Product Name: 14-3-3 β Rabbit Polyclonal Antibody**Catalog #: APRab06278**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	

Antigen Information


Gene Name	YWHAB
Alternative Names	YWHAB; 14-3-3 protein beta/alpha; Protein 1054; Protein kinase C inhibitor protein 1; KCIP-1
Gene ID	7529.0
SwissProt ID	P31946
Immunogen	The antiserum was produced against synthesized peptide derived from human 14-3-3 beta. AA range:41-90

Background


Research Area

Neuroscience


Image Data




Western blot analysis of lysates from HepG2 and Jurkat cells, using 14-3-3 beta Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using 14-3-3 β antibody.



Western Blot analysis of various cells using 14-3-3 β Polyclonal Antibody



Western Blot analysis of Jurkat cells using 14-3-3 β Polyclonal Antibody