
Product Name: ABHD12B Rabbit Polyclonal Antibody**Catalog #: APRab06434**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
Molecular Weight	40kDa

Antigen Information


Gene Name	ABHD12B
Alternative Names	ABHD12B; C14orf29; Abhydrolase domain-containing protein 12B
Gene ID	145447.0
SwissProt ID	Q7Z5M8
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD12B. AA range:281-330

Background

similarity:Belongs to the serine esterase family.,similarity:Belongs to the serine esterase family.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using ABHD12B Antibody. The picture on the right is blocked with the synthesized peptide.