

---

**Product Name: AK1 Rabbit Polyclonal Antibody****Catalog #: APRab06710**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	IHC,ICC/IF,ELISA
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

**Application**

**Dilution Ratio** IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000

**Molecular Weight**

**Antigen Information**

<b>Gene Name</b>	AK1
<b>Alternative Names</b>	AK1; Adenylate kinase isoenzyme 1; AK 1; ATP-AMP transphosphorylase 1; Myokinase
<b>Gene ID</b>	203.0
<b>SwissProt ID</b>	P00568
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human KAD1 . AA range:101-150

**Background**


adenylate kinase 1(AK1) Homo sapiens This gene encodes an adenylate kinase enzyme involved in energy metabolism and

homeostasis of cellular adenine nucleotide ratios in different intracellular compartments. This gene is highly expressed in skeletal muscle, brain and erythrocytes. Certain mutations in this gene resulting in a functionally inadequate enzyme are associated with a rare genetic disorder causing nonspherocytic hemolytic anemia. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015],catalytic activity:ATP + AMP = 2 ADP.,disease:Defects in AK1 are the cause of hemolytic anemia due to adenylate kinase deficiency [MIM:612631].,function:Catalyzes the reversible transfer of the terminal phosphate group between ATP and AMP. Small ubiquitous enzyme involved in energy metabolism and nucleotide synthesis that is essential for maintenance and cell growth.,online information:Adenylate kinase entry,polymorphism:This enzyme represents the most common of at least five alleles.,similarity:Belongs to the adenylate kinase family.,subunit:Monomer.,


## Research Area

Purine metabolism;

## Image Data



Immunofluorescence analysis of HepG2 cells, using KAD1 Antibody . The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using KAD1 Antibody . The picture on the right is blocked with the synthesized peptide.