
Product Name: ALS2CR7 Rabbit Polyclonal Antibody**Catalog #: APRab06814**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	46kDa

Antigen Information


Gene Name	CDK15 CDK15; ALS2CR7; PFTK2; Cyclin-dependent kinase 15; Amyotrophic lateral sclerosis 2
Alternative Names	chromosomal region candidate gene 7 protein; Cell division protein kinase 15; Serine/threonine-protein kinase ALS2CR7; Serine/threonine-protein kinase PFTAIR
Gene ID	65061.0
SwissProt ID	Q96Q40
Immunogen	The antiserum was produced against synthesized peptide derived from human ALS2CR7. AA range:261-310

Background

CDK15 (Cyclin-Dependent Kinase 15) is a Protein Coding gene. Diseases associated with CDK15 include amyotrophic lateral sclerosis 2, juvenile. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is CDK14. protein amino acid phosphorylation, phosphorus metabolic process, phosphate metabolic process, phosphorylation,

Research Area

Image Data



Western blot analysis of lysates from Jurkat and HUVEC cells, using AL2S7 Antibody. The lane on the right is blocked with the synthesized peptide.