

Summary

Production Name	ARMC6 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ARMC6
Alternative Names	ARMC6; Armadillo repeat-containing protein 6
Gene ID	93436.0
SwissProt ID	Q6NXE6.The antiserum was produced against synthesized peptide derived from human ARMC6. AA range:431-480

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:40000.Not yet tested in other applications.
Molecular Weight	54kDa

Product Name: ARMC6 Rabbit Polyclonal Antibody
Catalog #: APRab07147




Background


The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010],PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the ARMC6 family.,similarity:Contains 4 ARM repeats.,

Research Area

Image Data




Immunofluorescence analysis of A549 cells, using ARMC6 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HeLa and COLO cells, using ARMC6 Antibody. The lane on the right is blocked with the synthesized peptide.

Product Name: ARMC6 Rabbit Polyclonal Antibody
Catalog #: APRab07147



Western Blot analysis of various cells using ARMC6 Polyclonal Antibody

Note

For research use only.