
Product Name: Arrdc4 Rabbit Polyclonal Antibody**Catalog #: APRab07168**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
Molecular Weight	45kDa

Antigen Information


Gene Name	ARRDC4
Alternative Names	ARRDC4; Arrestin domain-containing protein 4
Gene ID	91947.0
SwissProt ID	Q8NCT1
Immunogen	The antiserum was produced against synthesized peptide derived from human ARRDC4. AA range:127-176

Background

similarity:Belongs to the arrestin family.,similarity:Belongs to the arrestin family.,

Research Area

Image Data



Western blot analysis of lysates from COLO cells, using ARR4 Antibody. The lane on the right is blocked with the synthesized peptide.