
Product Name: Ataxin-7L1 Rabbit Polyclonal Antibody**Catalog #: APRab07254**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000

Molecular Weight

Antigen Information


Gene Name	ATXN7L1
Alternative Names	ATXN7L1; ATXN7L4; KIAA1218; Ataxin-7-like protein 1; Ataxin-7-like protein 4
Gene ID	222255.0
SwissProt ID	Q9ULK2
Immunogen	The antiserum was produced against synthesized peptide derived from human ATXN7L1. AA range:441-490

Background

similarity:Contains 1 SCA7 domain.,similarity:Contains 1 SCA7 domain.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using ATXN7L1 Antibody. The picture on the right is blocked with the synthesized peptide.