
Product Name: Atg9a Rabbit Polyclonal Antibody**Catalog #: APRab07302**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:10000-1:20000
Molecular Weight	85kDa

Antigen Information


Gene Name	ATG9A APG9L1
Alternative Names	ATG9 autophagy related 9 homolog A (S. cerevisiae)
Gene ID	79065.0
SwissProt ID	Q7Z3C6
Immunogen	Synthetic Peptide of Atg9a

Background

function:Plays a role in autophagy.,similarity:Belongs to the ATG9 family.,function:Plays a role in autophagy.,similarity:Belongs to the ATG9 family.,

Research Area

Image Data



Western Blot analysis of 293T 3T3 cells using Atg9a Polyclonal Antibody diluted at 1:800. Secondary antibody was diluted at 1:20000