
Product Name: BTF3L4 Rabbit Polyclonal Antibody**Catalog #: APRab07687**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
Molecular Weight	17kDa

Antigen Information

Gene Name	BTF3L4
Alternative Names	BTF3L4; Transcription factor BTF3 homolog 4; Basic transcription factor 3-like 4
Gene ID	91408.0
SwissProt ID	Q96K17
Immunogen	The antiserum was produced against synthesized peptide derived from human BTF3L4. AA range:94-143


Background

similarity:Belongs to the NAC-beta family.,similarity:Contains 1 NAC-A/B (NAC-alpha/beta) domain.,similarity:Belongs to the

NAC-beta family.,similarity:Contains 1 NAC-A/B (NAC-alpha/beta) domain.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using BTF3L4 Antibody. The picture on the right is blocked with the synthesized peptide.