
Product Name: CaBP7 Rabbit Polyclonal Antibody**Catalog #: APRab07809**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
Molecular Weight	25kDa

Antigen Information


Gene Name	CABP7
Alternative Names	CABP7; CALN2; Calcium-binding protein 7; CaBP7; Calneuron II; Calneuron-2
Gene ID	164633.0
SwissProt ID	Q86V35
Immunogen	The antiserum was produced against synthesized peptide derived from human CABP7. AA range:1-50

Background

similarity:Contains 2 EF-hand domains.,similarity:Contains 2 EF-hand domains.,

Research Area

Image Data



Western Blot analysis of COLO cells using CaBP7 Polyclonal Antibody diluted at 1 : 1000