

Product Name: CLCC1 Rabbit Polyclonal Antibody
Catalog #: APRab08930



Summary

Production Name	CLCC1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	CLCC1
Alternative Names	CLCC1; KIAA0761; MCLC; Chloride channel CLIC-like protein 1; Mid-1-related chloride channel protein 1
Gene ID	23155.0
SwissProt ID	Q96S66.The antiserum was produced against synthesized peptide derived from human CLCC1. AA range:391-440

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:40000.Not yet tested in other applications.
-----------------------	---

Product Name: CLCC1 Rabbit Polyclonal Antibody
Catalog #: APRab08930




Molecular Weight 62kDa

Background


function: Seems to act as a chloride ion channel., similarity: Belongs to the chloride channel MCLC family., function: Seems to act as a chloride ion channel., similarity: Belongs to the chloride channel MCLC family.,

Research Area

Image Data




Immunofluorescence analysis of A549 cells, using CLCC1 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from COLO and HepG2 cells, using CLCC1 Antibody. The lane on the right is blocked with the synthesized peptide.


Product Name: CLCC1 Rabbit Polyclonal Antibody
Catalog #: APRab08930



Western blot analysis of the lysates from HT-29 cells using CLCC1 antibody.



Western Blot analysis of various cells using CLCC1 Polyclonal Antibody diluted at 1: 1000



Western Blot analysis of HepG2 cells using CLCC1 Polyclonal Antibody diluted at 1: 1000

Note

For research use only.