
Product Name: CY561 Rabbit Polyclonal Antibody**Catalog #: APRab09569**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	27kDa

Antigen Information


Gene Name	CYB561
Alternative Names	
Gene ID	1534.0
SwissProt ID	P49447
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110

Background

cofactor: Binds 2 heme groups non-covalently.,function: Secretory vesicle-specific electron transport protein.,similarity: Contains 1 cytochrome b561 domain.,cofactor: Binds 2 heme groups non-covalently.,function: Secretory vesicle-specific electron transport protein.,similarity: Contains 1 cytochrome b561 domain.,

Research Area

Image Data



Western blot analysis of lysates from U2OS cells, CY561 Rabbit Polyclonal Antibody was diluted at 1:1000, 4° over night